

WARRANTY
5
YEARS

CE

Made in France

200 kg

75 kg

140 mm

360/815 mm

2050 mm

997 mm

57°

33°



CONVERTISS

It is equipped with a translation raising backrest which aims to raise the patient in a secure way in order the best position while avoiding abdominal shear, and with a leg rest section with knee bending providing an optimized sitting comfort, reducing friction and preventing the patient from sliding down the bed.

The main advantage of the CONVERTISS model is its ability to easily transition, with affordable kits, from a standard medical bed (90 cm width - 200 cm length) to a more specialized model (up to 120 cm width - 220 cm length) to fit different morphologies without purchasing a new medical bed.

WATCH

the CONVERTISS product video on our YOUTUBE channel




YouTube
SCAN
to watch



Wooden bed panel available in 3 widths: 90 or 110 or 120 cm



Patient lifting pole in order to easily stand up with a stable point of support



2 enlargement kits available to obtain a width of 110 or 120 cm



Extension kit to obtain a length of 220 cm

REFERENCE
Convertiss: CONVERTISS-3F-PF

MEDICAL BEDS