

# CONVERTISS

## Bariatric XL

The CONVERTISS Bariatric XL provides 3 functions with a weight capacity of 200 kg as standard to handle bariatric patients.

It is equipped with a translation raising backrest which aims to raise the bariatric patient in a secure way in order the best position while avoiding abdominal shear, and with a leg rest section with knee bending providing an optimized sitting comfort, reducing friction and preventing the patient from sliding down.

Its main advantage is its ability to easily transition, with affordable kits, from a standard medical bed (90 cm width - 200 cm length) to a more specialized model (up to 120 cm width - 220 cm length) to fit different morphologies without purchasing a new medical bed.



REFERENCE CONVERTISS-3F-PF

DEMONSTRATION VIDEO

You Tube  
SCAN  
to watch



### KEY FEATURES



Wooden bed panel



Pair of folding side rails



Patient lifting pole in order to easily stand up with a stable point of support

### OPTIONS



Enlargement kits available to obtain a width of 110 or 120 cm

Extension kit to obtain a length of 220 cm

### TECHNICAL DATA

- A. Overall length: 206 cm
- B. Overall width: 90 cm
- C. Bed base length: 203 cm
- D. Bed base width: 84,5 cm
- E. Minimum height: 30,5 cm
- F. Maximum height: 80 cm

