

PRODUCT SHEET NAUSIFLOW 834-120-20 MATTRESS



Specifications:

Dimensions: 20 x 120 x 200 cm

20 cm air height

Top separable cover for easy maintenance -
CRIB 5 classified: antibacterial and antifungal
Class 1 medical device

Key features:

CPR (cardiopulmonary resuscitation): deflates quickly (less than 20 seconds with a patient) the portion of the mattress which supports the patient's bust

Top separable cover

Separable zipper

Independent cells in polyurethane on nylon frame

Quick fasteners for the disconnectable and removable cells

Transport mode