



300 kg



WARRANTY  
3 YEARS



Decubitus Level



SCAN ME  
to learn more

# NAUSIFLOW 100 MAXI system

## NAUSIFLOW 100 MAXI Pump

It manages bariatric patients up to 300 kg, operates in alternating/static modes, is designed for the care of stage 4 pressure ulcers and is used for prevention on very high-risk levels.

Its high-speed pump can manage a wide range of mattresses (air height up to 20 cm and width up to 140 cm) and it features a static operating mode to ensure treatment in palliative care.

## NAUSIFLOW 834-120 Air Mattress

It is adapted to the single sleeping width of 120 cm (20 x 120 x 200 cm - 20 cm full air height) for bariatric patients up to 300 kg.

It is equipped with 8-shaped cells which favours patient immersion by reducing contact pressure.



CPR (cardiopulmonary resuscitation): deflates quickly (less than 20 seconds with a patient) the portion of the mattress which supports the patient's bust



Transport mode



Timed care position



Bariatric program for patients above 150 kg

## REFERENCES

NAUSIFLOW 100 MAXI PUMP: **NA100MAXI-COMP-PM**  
NAUSIFLOW 834-120-20 MATTRESS: **NA834-MAT120-PM**  
NAUSIFLOW 100 MAXI SYSTEM 834-120: **NA100MAXI834120SYS**