

EASYLEV

mechanical opening legs



SCAN ME
to learn more

Its mechanical opening base legs (up to 94 cm) enables an ease of access to armchairs or wheelchairs besides the combined movement providing a smooth lifting process to the patient.

WATCH

the EASYLEV range product video
on our YOUTUBE channel



OPTIONS



Lithium Battery



Calf Strap

REFERENCES

with standard battery: EASYLEVECP-M-PF
with lithium battery: EASYLEVECP-M-LI-PF



150 kg



41 kg



690/1620 mm



113 mm



540/940 mm



Traction belt for a patented
combined movement in-
creasing patient comfort



Premium control box for
an enhanced use



Foam shin block for a
comfortable support



Mechanical opening legs
for better access

