



*Prior to each use, it is important to check the status of the device:  
seams, state of the straps and of the fabric.  
Any damage will either be repaired by NAUSICAA Medical, or will result in a disposal.*

## HEEL PADS



THIS MEDICAL DEVICE MUST BE  
IMPLEMENTED BY A TRAINED  
HEALTH PROFESSIONAL



### HEEL PADS SMALL MODEL


**HELP TO PREVENT  
PRESSURE SORES  
ON HEELS**

### HEEL PADS LARGE MODEL

**HELP TO PREVENT  
PRESSURE SORES  
ON HEELS**

#### Description :

The Heel Pads allow to to relieve pressure points and to wrap the foot in a furtrim, preventing rubbings on the heels and malleoli,



LMTAL-PF

### Important Recommendations :

In order to get a maximum efficiency of these sort of product, it is necessary :

- to choose the appropriate size for the patient
- to get the best possible adjustment to the patient

These products must not be in direct contact with a wounded skin.



## HELL PADS

COMBINATION OF POLYURETHAN FOAM  
AND JERSEY/COTTON POLYESTER FABRIC  
JERSEY 100% POLYESTER FABRIC  
100% POLYESTER FUR  
100% POLYESTER VELCRO STRAP

### INSTRUCTIONS FOR USE HELP TO PREVENT PRESSURE SORES ON HEELS

- Put the heel, fur side, right in the bottom of the heelpad.
- Close the heelpad on the foot through the velcro strap(s).

v 05.12/2020



THIS MEDICAL DEVICE MUST BE SET UP BY A PROFESSIONAL  
OR THE USER MUST BE EDUCATED TO SET IT UP.