

-  300 kg
-  42 kg
-  850/1350 mm
-  100 mm
-  690/1040 mm




SCAN ME  
to learn more

# WAYUP XL

electrical opening legs

The efficiency of its specific movement in one quick motion combined with its overpowering lifting actuator (10 000 N) allows the bariatric patient to reach a standing position while being transferred by activating the bust tonicity and supporting the lack of tonicity in the lower body.

**WATCH**  
the product video on  
our YOUTUBE channel



SCAN  
to watch






Premium control box for an enhanced use



Specially designed foam shin block




Electrical opening legs for ease of use



Removable platform to use the device as a walking aid

### OPTIONS



REFERENCES  
with standard battery: WAYUPXL-E-PF  
with lithium battery: WAYUPXL-E-LI-PF