

- 180 kg
- 40/42 kg
- 570/1780 mm
- 85 mm
- 630/1050 mm



SCAN ME  
to learn more

# FLYER PLUS

opening legs  
(mechanical or electrical)

It is specially designed for daily transfers due to its enhanced mobility (double ball bearings castors), mechanical or electrical opening (up to 105 cm), 8 000 N actuator (180 kg weight capacity) and our premium control box.

**WATCH**  
the FLYEVER/FLYER range product video  
on our YOUTUBE channel

VIDEO  
YouTube  
SCAN  
to watch



4-point spreader bar for an optimized weight distribution

Premium control box for an enhanced use

Electrical or mechanical opening legs for ease of use

75 mm castors with reduced height to pass under the lowest beds

### OPTIONS



REFERENCES MECHANICAL OPENING LEGS  
with standard battery: FLYERPLUS2-M-PF  
with lithium battery: FLYERPLUS2-M-LI-PF

REFERENCES ELECTRICAL OPENING LEGS  
with standard battery: FLYERPLUS2-E-PF  
with lithium battery: FLYERPLUS2-E-LI-PF