

EASYLEV XL

With Electrical Opening Legs

Made in France
CE
 WARRANTY 5 YEARS

780/1610 mm
 110 mm
 700/940 mm
 300 kg
 45 kg



KEY FEATURES



Traction belt for a combined movement increasing patient comfort



Electrical opening legs for ease of use



Removable platform to use the device as a walking aid

PRESENTATION

The EASYLEV is also available in XL version as a specialized model for bariatric patients up to 300 kg. Its ingenious patented combined movement between the telescopic mast and the traction sling allows the patient to reach a standing position thus delivering a high-quality care for bariatric patient.

OPTION



Lithium Battery for more efficiency:
 - up to 312 cycles
 - 2 years warranty

RECOMMENDED SLING

The Standard Sling (page 39)

Placed under the armpits, its full padding under the arms distributes traction on the patient thus avoiding painful pressure and supports the chest area.



REFERENCES

With standard battery	EASYLEV-XL-E-PF
With lithium battery	EASYLEV-XL-E-LI-PF