



NAUSI-FLY RANGE PRESENTATION



THOUGHT FOR OPTIMIZED COMFORT AND COST EFFECTIVENESS

With its great maneuverability, the NAUSI-FLY range has been conceived as an optimized mix between an economic budget and high-end specifications to perform secure and comfortable transfers with its smooth lifting movement.

It offers a simplicity of use and a versatility allowing to perform secure and comfortable transfers for all patients.

The NAUSI-FLY range is declined in 3 different models:

- Standard sized with mechanical opening or electrical opening (page 52)
- Compact sized with mechanical opening (page 53)

With the NAUSI-FLY range, a single caregiver can manage a wide range of patient handling routines in a secure environment while maintaining a correct posture to avoid injury.



LITHIUM BATTERY

As options, a lithium battery is available to enhance the use of the device and a mountable NAUSICAA weighing system (page 55) to determine the patient's weight.

NAUSICAA
Medical
NAUSI-FLY RANGE

SCAN
to watch



NAUSI-FLY

With Mechanical Opening Legs



NAUSI-FLY

With Electrical Opening Legs



**NAUSI-FLY
COMPACT**

With Mechanical Opening Legs

NAUSI-FLY RANGE KEY FEATURES

WARRANTY
5
YEARS



Equipped with advanced functional key features.

The NAUSI-FLY range is designed to have the best maneuverability to guarantee comfort to the patient and ease of use to the caregiver.

With their double large castors, these devices are fully adapted for home care and make it exceptionally easy to move around multiple floor types.

Moreover, they can lift a patient from the floor (page 49).



01

4-POINT WIDE SPREADER BAR for an optimized weight distribution for patient comfort
Better for:


- Large or tall patients
- Patient positioning at 90° for transfer to a chair
- Include head protection
- Lifting from floor



02

REMOTE AND CONTROL BOX

- SIMPLE FUNCTIONS or 4 FUNCTIONS (when electrical opening legs is chosen) with battery display level on the remote control



03

OPENING LEGS:

- Mechanical
- Electrical



04

CASTORS:

- 4 Double castors and 2 double braked castors 100 mm
- Frame height : 113 mm

