

EASYLEV XL

electrical opening legs



The specificity of its doubled movement combined with its overpowering lifting actuator (10 000 N) allows caregivers to provide high quality care during transfers as it supports a very low tonicity of the bariatric patient's bust during the standing-up process.

WATCH

the product video on
our YOUTUBE channel



- 300 kg
- 45 kg
- 730/1650 mm
- 113 mm
- 690/1040 mm

OPTIONS



Lithium Battery



Calf Strap




Traction belt for a combined movement increasing patient comfort



Premium control box for an enhanced use



Electrical opening legs for ease of use



Removable platform to use the device as a walking aid



REFERENCES

with standard battery: EASYLEV-XL-E-PF
with lithium battery: EASYLEV-XL-E-LI-PF