

NAUSI-FLY COMPACT

With Mechanical Opening Legs



CE
WARRANTY 5 YEARS

580/1600 mm

110 mm

570/880 mm

160 kg

31 kg



KEY FEATURES



4-point spreader bar for an optimized weight distribution



Large handlebar for ease of handling



4 large doubles castors with double ball bearing for improved mobility

PRESENTATION

Specially designed for home care use where space is limited, its reduced dimensions and its 160 kg weight capacity (6 000 N actuator) make it a remarkably efficient patient lift.

OPTION



Lithium Battery for more efficiency:
- up to 84 cycles
- 2 years warranty

REFERENCES

With standard battery	NAFLY-160-CPT-PF
With lithium battery	NAFLY-160-CP-LI-PF

