

# FLYER PLUS

With Mechanical or Electrical Opening Legs



## KEY FEATURES



4-point spreader bar for an optimized weight distribution



Premium control box for an enhanced use



75 mm castors with reduced height to pass under the lowest beds

## PRESENTATION

Our high-end device is specially conceived for daily transfers due to its enhanced mobility (double ball bearings castors), its mechanical (up to 109 cm) or electrical (up to 114 cm) opening legs and its 8 000 N actuator (180 kg weight capacity).

## OPTIONS



Lithium Battery for more efficiency:  
- up to 84 cycles  
- 2 years warranty



A mountable NAUSICAA weighing system to determine the patients weight (page 55)

## REFERENCES

Mechanical Opening	With standard battery	FLYERPLUS2-M-PF
	With lithium battery	FLYERPLUS2-M-LI-PF
Electrical Opening	With standard battery	FLYER+2-E2-PF
	With lithium battery	FLYER+2-E2-LI-PF

