

# NAUSIFLOW 100 AUTO QUATTRO System



150 kg

WARRANTY  
**3**  
YEARS



**4**  
Decubitus  
Level



## NAUSIFLOW 100 AUTO QUATTRO Pump

It operates in alternating/static modes up to 150 kg, is designed for the care of stage 4 pressure ulcers and is used for prevention on very high-risk levels.

It automatically adjusts to any air level, is equipped with a key lock, a timed care position, a classic dynamic mode, a secured static mode in low continuous pressure, a sitting position and 3 comfort programs.



### KEY FEATURES

- 01** 4 different modes of air flow
- 02** 4 different programs of pressure level
- 03** Timed care position
- 04** Panel lock

## NAUSIFLOW 100-512 - NAUSIFLOW 100-834 Air Mattresses

All information on the NAUSIFLOW 100-512 Air Mattress can be found on page 86. The NAUSIFLOW 100-834 Air Mattress is equipped with 8-shaped cells which enable patient immersion by reducing contact pressure and it allows multiple configurations due to its dischargeable cells with plug.

Its compartmented configuration with a double layer of air and its adjusted size (20 x 85 x 200 cm – 20 cm full air height) fit any type of 90 cm width bed (box spring from 85 to 90 cm).



Independent cells in polyurethane on nylon frame



Transport mode



CPR

REFERENCES	
NAUSIFLOW AUTO 100 QUATTRO Pump	NA100QTO-COMP-PM
NAUSIFLOW 100-512 Mattress	NA512-MAT85-PM
NAUSIFLOW 100-834 Mattress	NA834-MAT85-PM
NAUSIFLOW 100 AUTO System 512	NA100QTO512-SYS
NAUSIFLOW 100 AUTO System 834	NA100QTO834-SYS