

# NAUSIFLOW 100 AUTO QUATTRO

system



SCAN ME  
to learn more

## NAUSIFLOW 100 AUTO QUATTRO Pump

It operates in alternating/static modes up to 150 kg, is designed for the care of stage 4 pressure ulcers and is used for prevention on very high-risk levels.

It automatically adjusts to any air level, is equipped with a key lock, a timed care position, a classic dynamic mode, a secured static mode in low continuous pressure, a sitting position and 3 comfort programs.

## NAUSIFLOW 100-512 & NAUSIFLOW 100-834 Air Mattress

All information on the NAUSIFLOW 100-512 Air Mattress can be found on page 75.

The NAUSIFLOW 100-834 Air Mattress is equipped with 8-shaped cells which favours patient immersion by reducing contact pressure and it enables multiple configurations due to its dischargeable cells with plug.

Its compartmented configuration with a double layer of air and its adjusted size (20 x 85 x 200 cm – 20 cm full air height) fit any type of 90 cm width bed (box spring from 85 to 90 cm).



150 kg



WARRANTY  
3  
YEARS



Decubitus  
Level



PRESSURE CARE SYSTEMS

### REFERENCES

NAUSIFLOW AUTO 100 QUATTRO PUMP: **NA100QTO-COMP-PM**  
 NAUSIFLOW 100-512 MATTRESS: **NA512-MAT85-PM**  
 NAUSIFLOW 100-834 MATTRESS: **NA834-MAT85-PM**  
 NAUSIFLOW 100 AUTO SYSTEM 512: **NA100QTO512-SYS**  
 NAUSIFLOW 100 AUTO SYSTEM 834: **NA100QTO834-SYS**

CPR (cardiopulmonary resuscitation): deflates quickly (less than 20 seconds with a patient) the portion of the mattress which supports the patient's bust

Independent cells in polyurethane on nylon frame

4 different modes of air flow

4 different programs of pressure level