

ACTIV'UP

With Mechanical Opening Legs

Made in France
CE
 WARRANTY 5 YEARS

150 kg
 29 kg



The ACTIV'UP is a versatile device for active transfers as it is:



1. A lifting aid handlebar



2. A secure seated or standing transfer system



3. A walking aid

PRESENTATION

The ACTIV'UP is used in rehabilitation to stimulate the patient during a transfer.

It encourages patients to participate in their daily transfers, allowing them to improve their blood circulation, breathing, transit and muscle tone.

REFERENCE

ACTIVUP-PF

TECHNICAL DATA

- A. Total length: 100 cm
- B. Minimum lifting height: 75 cm
- C. Maximum lifting height: 89 cm (with extension)
- D. Total height: 145 cm
- E. Height of chassis: 11 cm
- F. Minimum base width: 57 cm
- G. Maximum base width: 88 cm
- H. Turning diameter: 112 cm

