

The ACTIV'UP is a versatile device for active transfers as it is:



1. A lifting aid handlebar




2. A secure seated or standing transfer system



3. A walking aid

# ACTIV'UP SLING & EXTENSION



1. The ACTIV'UP sling is specially designed to be used with the ACTIV'UP to help the patient during the lifting process and can serve as a back panel for sitting positions.

2. The extension raises the mast height of 14 cm for taller patients when used as a walking aid.

REFERENCES	ACTIV'UP SLING	ACT-SGL-PF
	ACTIV'UP EXTENTION	ACT-RAL