

ACTIV'UP

With Mechanical Opening Legs

The ACTIV'UP is used in rehabilitation to stimulate the patient during a transfer.

It encourages patients to participate in their daily transfers, allowing them to improve their blood circulation, breathing, transit and muscle tone.



REFERENCE

ACTIVUP-PF

DEMONSTRATION VIDEO

YouTube
SCAN
to watch



KEY FEATURES



Removable platform to use the device as a walking aid



Ergonomic handlebar to help the patient stand



Fabric seating for the patient's comfort

TECHNICAL DATA

- A. Total length: 100 cm
- B. Minimum lifting height: 75 cm
- C. Maximum lifting height: 89 cm (with extension)
- D. Total height: 145 cm
- E. Height of chassis: 11 cm
- F. Minimum base width: 57 cm
- G. Maximum base width: 88 cm
- H. Turning diameter: 112 cm

