

*It is important to check the condition of the device: seams, state of the fabric.
Any deterioration will result in either Nausicaa Medical repair or scrapping.*

SLIDING GLOVES



**THIS MEDICAL DEVICE MUST BE
IMPLEMENTED OR TRAINED BY A
HEALTHCARE PROFESSIONAL**



MDG-PF

SLIDING GLOVES

HOSPITAL RANGE
100% polyamide warp and weft
1-sided silicone coating

DESCRIPTION

- Manual patient transfers are an integral part of the working day for caregivers.
- For the handling of a patient, the Glide Sheet is a device that facilitates the lateral and / or vertical translation as well as the raising of a patient lying in a bed, thus limiting back pain and musculoskeletal disorders.
- Sliding effect facilitates the work of caregivers by reducing friction due to patient weight during handling.
- It consists of a synthetic polyester fabric tube (polyamide with silicone-treated polyurethane coating), the two faces of fabric in sliding contact one on the other.

v 03.12/2020



THIS MEDICAL DEVICE MUST BE IMPLEMENTED OR MAKE
THE SUBJECT OF TRAINING BY A HEALTHCARE PROFESSIONAL