

WAYUP 5

With Fixed Base

Made in France
CE
 WARRANTY 5 YEARS

800/1330 mm
 65 mm
 240/565 mm
 180 kg
 29 kg



KEY FEATURES



Innovative expandable lifting arm



Ergonomic knee block to protect vulnerable knees



Compact base adapted to toilet shape

PRESENTATION

Totally different than other devices on the market, our best-seller is specially optimized for confined spaces with its extreme compactness (56 cm base width) and its remarkably low toilet shape base (6,5 cm overall height). Its weight capacity has been increased to 180 kg.

OPTION



Lithium Battery for more efficiency:
 - up to 312 cycles
 - 2 years warranty

RECOMMENDED SLING

The Specific Sling (page 38)

Positioned on the lumbar region in order to trigger the patient's natural raising scheme, its rigid foam padding at the back-part level reduces painful pressure and its belt allows transfers in a sitting position.



REFERENCES

With standard battery	WAYUP5-PF
With lithium battery	WAYUP5-LI-PF

