

EASYLEV 5

With Fixed Base

Made in France
CE
 WARRANTY 5 YEARS

720/1620 mm
 65 mm
 240/565 mm
 180 kg
 38 kg



KEY FEATURES



Traction belt for a combined movement increasing patient comfort



Innovating extendable lifting arms



Compact base adapted to toilet shape

PRESENTATION

In addition to its combined movement which is specifically adapted for patients with very low tonicity in the chest area, it is equipped with a special shape to be adapted to a standard toilet. It is also perfectly adapted to confined spaces with its remarkable low base (6,5 cm overall height) and its extraordinary compactness (56 cm base width).

OPTION



Lithium Battery for more efficiency:
 - up to 312 cycles
 - 2 years warranty

RECOMMENDED SLING

The Standard Sling (page 39)

Placed under the armpits, its full padding under the arms distributes traction on the patient thus avoiding painful pressure and supports the chest area.



REFERENCES

With standard battery	EASYLEV5-PF
With lithium battery	EASYLEV5-LI-PF