



SCAN ME
to view all
devices

STAND-UP LIFTS

Presentation

All our stand-up lifts are compliant to European Standard NF EN ISO 10535:2007.

Stand-up lifts make caregivers' work easier

Stand-up lifts are designed to enable patients to be raised up from a bed, a comfort chair, a toilet seat or a wheelchair for safe transfers.

They allow a single caregiver to lift a patient with minimum effort in order to devote more attention for encouraging the patient's active participation in the transfer routine.

The process is fast, safe, and efficient

Transfer with a stand-up lift is qualified as a "semi-active transfer" because the patient needs to participate during the movement and maintain some stimulus (in at least one of the two legs).

Encouraging patient mobility and promoting dignity, stand-up lifts are ideal for patients whose mobility is reduced and have difficulty getting up into a standing position.

Patients in home or in care environments gain significant physical and psychological benefits when their mobility is preserved.

What you need to know before choosing

Selecting the right device is essential for maintaining mobility as the use of adapted equipment reduces the risk of injury, both for the patient and for the caregiver.

5
YEARS

Warranty

CE

Certified



Quality



Lithium Battery
Available

