



EASYLEV RANGE

A HIGH-END LINE FOR THE BEST QUALITY CARE

Characterized by its ingenious patented double movement, the EASYLEV range is a high-end line as it allows caregivers provide high quality care with the best comfort in all environments.

Declined in 3 different models (fixed base or mechanical or electrical opening), it is also available in a bariatric version which includes the features of the EASYLEV 2 (page 24).

Extremely comfortable with its smoothly motion.

Due to its combined movement of the telescopic mast and of the traction sling, the EASYLEV range is providing a high level of comfort to the patient by supporting a very low tonicity of the patient's chest area during all the standing-up process.



As an option, a lithium battery (page 25) is available to enhance the use of the device.



NAUSICAA Medical EASYLEV RANGE



SCAN
to watch



WARRANTY

5 YEARS



EASYLEV 5 Fixed Base



EASYLEV 5 With Mechanical Opening Legs



EASYLEV 5 PLUS With Electrical Opening Legs



EASYLEV XL With Electrical Opening Legs



STAND-UP LIFTS