

KIT - WAYUP 4 TO WAYUP 5



Castors Kit



Actuator Kit



Lifting Arms Kit



Ergonomic Knee Block Kit

KIT - EASYLEV 2 TO EASYLEV 5



Castors Kit



Actuator Kit



Lifting Arms Kit



Ergonomic Knee Block Kit

WAYUP & EASYLEV KITS

As the new WAYUP and EASYLEV have arrived, NAUSICAA Medical offers you the possibility to upgrade your former devices with a kit containing the following parts :

| REFERENCES | | |
|------------|--------------------------|-------------|
| WAYUP | Castors Kit | KT-RB-VERT4 |
| | Actuator Kit | KT-VRN-WP5 |
| | Lifting Arms Kit | WP5-BRS-PF |
| | Ergonomic Knee Block Kit | KT-MAG |

| REFERENCES | | |
|------------|--------------------------|-------------|
| EASYLEV | Castors Kit | KT-RB-VERT4 |
| | Actuator Kit | KT-VRN-ESL5 |
| | Lifting Arms Kit | ESL5-BRS-PF |
| | Ergonomic Knee Block Kit | KT-MAG |

LEAD BATTERY

LITHIUM BATTERY

Autonomy on a 11 cm lifting stroke actuator (WAYUP, BLUE WAYUP)

130 cycles

312 cycles

Time to recharge up to 85%

2h30

1h10

Lifespan (number of complete charges)

Up to 300

Minimum of 2000

Lifespan (number of micro charges)

Up to 240

Up to more than 8000

Weight of the battery pack

4 kg

1,8 kg

Warranty

None

2 years

Average lifetime

6 months to 2 years

2 years to 8 years

LITHIUM BATTERY OPTION

More efficient at all levels and also lighter in handling, lithium batteries are perfectly suited for intensive use conditions.

Lithium batteries support very well partial charging which reflects more the reality of use as you can recharge them as soon as the device is no longer used, even for a few minutes.

