

# STAND-UP LIFTS Comparative



Weight Capacity  
(in kg)



Total Weight  
(in kg)



Minimum and  
Maximum Lifting  
Height (in cm)



Frame Height  
(in cm)



Internal and  
External Width  
Legs (in cm)



Minimum and  
Maximum Base  
Width (in cm)

	Weight Capacity (in kg)	Total Weight (in kg)	Minimum and Maximum Lifting Height (in cm)	Frame Height (in cm)	Internal and External Width Legs (in cm)	Minimum and Maximum Base Width (in cm)
WAYUP 5 Fixed Base (page 14)	180	29	80 / 133	6,5	24 / 56,5	-
WAYUP 5 Mechanical Opening Legs (page 15)	180	36	82 / 135	11	-	54 / 90
WAYUP 5 Electrical Opening Legs (page 15)	180	38	82 / 135	11	-	54 / 100
WAYUP XL Electrical Opening Legs (page 16)	300	45	92 / 131	11	-	70 / 94
EASYLEV 5 Fixed Base (pages 20)	180	38	72 / 162	6,5	24 / 56,5	-
EASYLEV 5 Mechanical Opening Legs (page 21)	180	42	75 / 165	11	-	54 / 90
EASYLEV 5 Electrical Opening Legs (page 21)	180	44	75 / 165	11	-	54 / 100
EASYLEV XL Electrical Opening Legs (page 22)	300	45	78 / 161	11	-	70 / 94
BLUE WAYUP Fixed Base (page 26)	150	29	79 / 123	8,5	29,5 / 59	-
BLUE WAYUP Mechanical Opening Legs (page 27)	150	33	84 / 128	11	-	57 / 88
BLUE WAYUP Electrical Opening Legs (page 27)	150	38,5	84 / 128	11	-	57 / 90
BLUE EASYLEV Fixed Base (page 32)	150	38	66 / 156	8,5	29,5 / 59	-
BLUE EASYLEV Mechanical Opening Legs (page 33)	150	42	70 / 160	11	-	57 / 88
BLUE EASYLEV Electrical Opening Legs (page 33)	150	44	70 / 160	11	-	57 / 90