


# CONVERTISS

Folding Premium



## KEY FEATURES



Quick electric folding and unfolding function



Patient lifting pole in order to easily stand up with a stable point of support



Pair of folding side rails for more safety

## PRESENTATION

The CONVERTISS Folding Premium offers the same functions as the CONVERTISS Standard Plus with the addition of a folding/unfolding electrical system for easy transport and installation along with an enhanced weight capacity of 180 kg.

It is a home care electric folding bed which has been specially designed to be installed by just one person; as simple as pressing a button on the remote control, the installer can unfold the bed very quickly.

Head and foot sections can be adjusted for the patient's comfort and the height can be adjusted from 70 cm to 30 cm from the floor.

As an option, in addition to a full-length side rails and a bed lever to enhance the use of the bed, an enlargement kit is available (see page 82) to extend the sleeping width up to 120 cm.

REFERENCE

NAUSIBED-6F