

CONVERTISS

Standard Plus

The CONVERTISS Standard Plus, whose maximum weight capacity is 135 kg, is equipped with 4 actuators allowing 5 different functions: variable height, raising backrest section, adjustable leg rest section, chair position and Trendelenburg position.

Easy to use, its space saving storage facility is the main interest to Nursing, Care Home and Community environments.

The Trendelenburg position (and reverse Trendelenburg) is made to optimize venous return to the heart and prevents pooling of blood in the lower extremities.

As an option, an enlargement kit (up to 120 cm), a full-length side rails and a bed lever are available to enhance the use of the bed (see page 73).



REFERENCE NAUSIBED-5F

NEW Product

DEMONSTRATION VIDEO

YouTube
SCAN
to watch



KEY FEATURES



Transport kit



Patient lifting pole in order to easily stand up with a stable point of support



Pair of folding side rails for more safety

TECHNICAL DATA

- A. Overall length: 215 cm
- B. Overall width: 99 cm
- C. Bed base length: 203 cm
- D. Bed base width: 83,5 cm
- E. Minimum height: 24 cm
- F. Maximum height: 65 cm
- G. Folded bed width: 45 cm
- H. Folded bed length: 99 cm
- I. Folded bed height: 135 cm

