


MEDICAL BEDS

Comparative



CONVERTISS
Standard Plus
page 70

CONVERTISS
Folding Premium
page 71

CONVERTISS
Bariatric XL
page 72

| | Capacity (in kg) | Lowest/Highest Height Position (in cm) | Working Length (in cm) | Working Width (in cm) | Extension Kit to 120 cm | Backrest | Knee Break | Side Rails | Full Length Side Rails | Patient Lifting Pole | Bed Weight (in kg) |
|------------------------------------|------------------|--|------------------------|-----------------------|-------------------------|------------|------------|------------|------------------------|----------------------|--------------------|
| CONVERTISS Standard Plus page 70 | 135 | 24/65 | 215 | 99 | Option | Electrical | Electrical | Yes | Option | Yes | 69 |
| CONVERTISS Folding Premium page 71 | 180 | 31/70 | 215 | 90 | Option | Electrical | Electrical | Yes | Option | Yes | 79 |
| CONVERTISS Bariatric XL page 72 | 200 | 30/80 | 206 | 90 | No | Electrical | Electrical | Yes | No | Yes | 75 |