

WAYUP 5

With Mechanical or Electrical Opening Legs



KEY FEATURES



Innovating extendable lifting arms



Ergonomic knee block to protect vulnerable knees



Removable low footrest for tailored uses



Mechanical or electrical opening legs for better access

PRESENTATION

This mechanical or electrical opening legs version makes it easy to fit around armchairs or wheelchairs. Both models are perfectly adapted to efficiently perform a change of the patient's lower body and can be used as a walking aid by removing the foot platform.

- 2 opening base legs versions:
- Mechanical opening legs
 - Electrical opening legs

OPTION



Lithium Battery for more efficiency:
 - up to 312 cycles
 - 2 years warranty

REFERENCES

Mechanical Opening	With standard battery	WAYUP5ECP-M-PF
	With lithium battery	WAYUP5ECP-M-LI-PF
Electrical Opening	With standard battery	WP5ECP-E2-PF
	With lithium battery	WP5ECP-E2-LI-PF

RECOMMENDED SLING

The Specific Sling (page 38)



Positioned on the lumbar region in order to trigger the patient's natural raising scheme, its rigid foam padding at the back-part level reduces painful pressure and its belt allows transfers in a sitting position.