

EASYLEV 5

With Mechanical or Electrical Opening Legs

Made in France
CE
 WARRANTY 5 YEARS

750/1650 mm
 110 mm
 M: 540/900 mm
 E: 540/1000 mm
 180 kg
 42/44 kg



KEY FEATURES



Innovating extendable lifting arms



Ergonomic knee block to protect vulnerable knees



Removable low footrest for tailored uses



Mechanical or electrical opening legs for better access

PRESENTATION

This mechanical or electrical opening legs version makes it easy to fit around armchairs or wheelchairs. Specially designed for patients who have a very low tonicity in the chest area, they are both completely dedicated to a smooth comfortable standing-up process due to their specific combined movement.

- 2 opening base legs versions:
- Mechanical opening legs
 - Electrical opening legs

OPTION



Lithium Battery for more efficiency:
 - up to 312 cycles
 - 2 years warranty

REFERENCES

Mechanical Opening	With standard battery	EASYLEV5ECP-M-PF
	With lithium battery	EASYLEV5ECP-M-LI-PF
Electrical Opening	With standard battery	ESL5ECP-E2-PF
	With lithium battery	ESL5ECP-E2-LI-PF

RECOMMENDED SLING

The Standard Sling (page 39)

Placed under the armpits, its full padding under the arms distributes traction on the patient thus avoiding painful pressure and supports the chest area.

