

BLUE WAYUP

With Mechanical or Electrical Opening Legs



KEY FEATURES



Fixed handles for a tight grip



Foam shin block for a comfortable support



Mechanical or electrical opening legs for better access



Removable platform to use the device as a walking aid

PRESENTATION

This mechanical or electrical opening legs version is specially designed to provide an optimal ease of access to armchairs or wheelchairs. They are perfectly adapted to efficiently perform a change of the patient's lower body and be easily transformed into a walking aid by removing the foot platform.

- 2 opening base legs versions:
- Mechanical opening legs
 - Electrical opening legs

OPTION



Lithium Battery for more efficiency:
 - up to 312 cycles
 - 2 years warranty

REFERENCES

Mechanical Opening	With standard battery	WAYUPBLUECP-PF
	With lithium battery	WAYUPBLUECP-LI-PF
Electrical Opening	With standard battery	WAYUPBLEP+-E-PF
	With lithium battery	WAYUPBLEP+-E-LI-PF

RECOMMENDED SLING

The Eco Sling (page 38)

Placed under the armpits or on the lumbar region, its padding under the arms distributes traction on the patient and avoids painful pressure during a transfer.

