

BLUE EASYLEV

With Mechanical or Electrical Opening Legs

Made in France
CE
WARRANTY 5 YEARS

700/1600 mm 110 mm M: 570/880 mm 150 kg 42/44 kg
 E: 570/900 mm



KEY FEATURES



Traction belt for a combined movement increasing patient comfort



Continuous handlebar for ease of handling and better maneuverability during the complete transfer



Mechanical or electrical opening legs for ease of use




Removable platform to use the device as a walking aid



PRESENTATION

This mechanical or electrical opening legs version is specially designed to provide an optimal ease of access to armchairs or wheelchairs. Both models are completely dedicated to a smooth comfortable standing-up process due to their specific combined movement. By removing the foot platform, they can be easily transformed into a walking aid to encourage the patient to walk on its own.

OPTION

 Lithium Battery for more efficiency:
 - up to 312 cycles
 - 2 years warranty

REFERENCES

Mechanical Opening	With standard battery	EASYLEVBLUECP-PF
	With lithium battery	EASYLEVBLUECP-LI-PF
Electrical Opening	With standard battery	ESLB2EP+-E2-PF
	With lithium battery	ESLB2EP+-E2-LI-PF

RECOMMENDED SLING



The Eco Sling (page 38)

Placed under the armpits or on the lumbar region, its padding under the arms distributes traction on the patient and avoids painful pressure during a transfer.