

# BLUE EASYLEV

With Fixed Base

Made in France  
**CE**  
 WARRANTY 5 YEARS

660/1560 mm  
 85 mm  
 295/590 mm  
 150 kg  
 38 kg



## KEY FEATURES



Traction belt for a combined movement increasing patient comfort



Continuous handlebar for ease of handling and better maneuverability during the complete transfer



Foam shin block for a comfortable support



Ball bearing castors for a better maneuverability

## PRESENTATION

It has been designed for patients with very low tonicity in the chest area. Easy to handle with its oversized handlebar, it is made for confined spaces with its extremely compact base (59 cm base width). It allows easy access for the caregiver to perform a change during the daily toilet.

## OPTION



Lithium Battery for more efficiency:  
 - up to 312 cycles  
 - 2 years warranty

## RECOMMENDED SLING

The Eco Sling (page 38)

Placed under the armpits or on the lumbar region, its padding under the arms distributes traction on the patient and avoids painful pressure during a transfer.



### REFERENCES

With standard battery	EASYLEVBLUE-PF
With lithium battery	EASYLEVBLUE-LI-PF