

# NAUSI-FLY

With Mechanical or Electrical Opening Legs



CE 5 YEARS WARRANTY

600/1750 mm

110 mm

M: 680/980 mm  
E: 680/1080 mm

180 kg

38/40 kg



## KEY FEATURES



4-point spreader bar for an optimized weight distribution



8 000 N actuator offering a compromise between maximum capacity (180 kg) and lifting speed



Electrical opening legs for ease of use



4 large doubles castors with double ball bearing for improved mobility

## PRESENTATION

These mechanical and electrical opening legs versions allow enhanced maneuverability for daily transfers with their double ball bearing castors. In addition to the 8 000 N actuator (180 kg weight capacity), the 4-points spreader bar is perfectly adapted to provide a maximum comfort and security to the patient.

## OPTIONS



Lithium Battery for more efficiency:  
- up to 84 cycles  
- 2 years warranty



A mountable NAUSICAA weighing system to determine the patients weight (page 55)

## REFERENCES

Mechanical Opening	With standard battery	NAFLY2-180-PF
	With lithium battery	NAFLY2-180-LI-PF
Electrical Opening	With standard battery	NAFLY+2-PF
	With lithium battery	NAFLY+2-LI-PF