

# BLUE FLYER

With Mechanical or Electrical Opening Legs

Made in France  
**CE**  
 WARRANTY 5 YEARS

650/1800 mm  
 110 mm  
 M: 680/980 mm  
 E: 680/1080 mm  
 180 kg  
 38/40 kg



## KEY FEATURES



360-degree rotation of the spreader bar for a simplified positioning



8 000 N actuator offering a compromise between maximum capacity (180 kg) and lifting speed



Electrical opening legs for ease of use



Reduced height of the front castors to pass under the lowest beds

## PRESENTATION

These mechanical and electrical opening legs versions are extremely easy to use and specially designed to perform secure and comfortable transfers for most patients (180 kg weight capacity with an 8 000 N actuator).

## OPTIONS



Lithium Battery for more efficiency:  
 - up to 84 cycles  
 - 2 years warranty



A mountable NAUSICAA weighing system to determine the patients weight (page 55)

## REFERENCES

		REFERENCES
Mechanical Opening	With standard battery	FLYER-BLUE-PF
	With lithium battery	FLYER-BLUE-LI-PF
Electrical Opening	With standard battery	FLYERBLUE+-E-PF
	With lithium battery	FLYERBLUE+-E-LI-PF