


# TECHNICAL DATA BLUE WAYUP RANGE

## BLUE WAYUP Fixed Base




### DIMENSIONS

- A. Total length: 94 cm
- B. Minimum lifting height: 79 cm
- C. Maximum lifting height: 123 cm
- D. Total height: 95 cm
- E. Frame height: 8,5 cm
- F. Minimum base width: 29,5 cm
- G. Maximum base width: 59 cm
- H. Turning diameter: 90 cm




## BLUE WAYUP With Mechanical Opening Legs



### DIMENSIONS

- A. Total length: 100 cm
- B. Minimum lifting height: 84 cm
- C. Maximum lifting height: 128 cm
- D. Total height: 100 cm
- E. Frame height: 11 cm
- F. Minimum base width: 57 cm
- G. Maximum base width: 88 cm
- H. Turning diameter: 115 cm



## BLUE WAYUP With Electrical Opening Legs



### DIMENSIONS

- A. Total length: 100 cm
- B. Minimum lifting height: 84 cm
- C. Maximum lifting height: 128 cm
- D. Total height: 100 cm
- E. Frame height: 11 cm
- F. Minimum base width: 57 cm
- G. Maximum base width: 90 cm
- H. Turning diameter: 115 cm

