


TECHNICAL DATA NAUSI-FLY RANGE

NAUSI-FLY With Mechanical Opening Legs



DIMENSIONS

- A. Total length: 119 cm
- B. Minimum lifting height: 60 cm
- C. Maximum lifting height: 175 cm
- D. Total height: 138 cm
- E. Frame height: 11 cm
- F. Minimum base width: 68 cm
- G. Maximum base width: 98 cm
- H. Turning diameter: 135 cm




NAUSI-FLY With Electrical Opening Legs



DIMENSIONS

- A. Total length: 119 cm
- B. Minimum lifting height: 60 cm
- C. Maximum lifting height: 175 cm
- D. Total height: 138 cm
- E. Frame height: 11 cm
- F. Minimum base width: 68 cm
- G. Maximum base width: 108 cm
- H. Turning diameter: 135 cm



NAUSI-FLY COMPACT With Mechanical Opening Legs



DIMENSIONS

- A. Total length: 100 cm
- B. Minimum lifting height: 58 cm
- C. Maximum lifting height: 160 cm
- D. Total height: 126 cm
- E. Frame height: 11 cm
- F. Minimum base width: 57 cm
- G. Maximum base width: 88 cm
- H. Turning diameter: 115 cm

